

Isolation, characterization and stabilization of urinary Extracellular vesicles (EVs) and EV RNA

Amit Arora^{1*}, Jaclyn Chabot², Gagandeep Narula¹, Evgueni Doukhanine¹, Anthoula Lazaris², Peter Metrakos², Vanessa Vankerckhoven³, Rafal Iwaszow¹

¹ DNA Genotek, 3000-500 Palladium Drive, Ottawa, ON Canada K2V 1C2

² Research Institute-McGill University Health Centre, Royal Victoria Hospital, 1001 Decarie Blvd, Montreal, Quebec, Canada H4A 3J1

³ Novosanis NV, Bijkhoevelaan 32c, 2110 Wijnegem, Belgium

Email: amit.arora@dnagenotek.com



EV RNA content in different urine fractions.

Urine fractions studied:

- First morning first void (FMFV)
- First morning first midstream (FMMS)
- Mid-day/ random first void (MDFV)

Extraction techniques used:

- Qiagen ExoRNeasy Maxi kit


RNA quantification:

- RT-qPCR assay on β -actin (ACTB)



Colli-Pee device (Novosanis)


Figure 1



DNA Genotek's proprietary urinary analyte stabilizer (UAS) preserves detection of EV RNA in urine samples for 7 days at RT.


- Mid-day/ random First void urine samples with and without UAS (UAS:urine ratio = 0.4:1).
- EV RNA extraction: Qiagen ExoRNeasy Maxi kit
- RNA Quantification: RT-qPCR assay on β -actin (ACTB)

Figure 2



Total EV concentration appears to be both donor and extraction method dependent.

Figure 3



2UC- Ultra-centrifugation; Exo- Qiagen Exo Easy kit; miR- Qiagen miRCURY exosome isolation kit

Evidence for urinary exosome concentration to be both donor and extraction method dependent.

Figure 4

