

Data Sheet

Colli-Pee® containing UCM preservative, FV-5000 series



Product Description

Product Name: Colli-Pee® containing UCM preservative, FV-5000 series

Reference Number: N00176

Components: Housing, floater, prefilled collector tube, tube cap and IFU

25 or 100 Order Quantity:

Product Characteristics

Material: Polypropylene (PP) and High Density Polyethylene (HDPE) Packaging: See-through pouch constructed of medical grade paper/

multiply PET/PP-plastic laminate

Tube Volume: 30 mL

Collected Volume: 18-22 mL 7 mL

Preservative Volume:

Collector Tube - Tube Cap Interface: Screw thread Housing - Collector Tube Interface: Screw thread

Shelf-life: 2 years Non-sterile Sterility:

