


# Data Sheet

Colli-Pee® Small Volumes containing UCM




## Product Description

Product Name:	Colli-Pee® Small Volumes containing UCM
Reference Number:	N00327
Components:	Housing, floater, prefilled collector tube, tube cap and IFU
Order Quantity:	25 or 100

## Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of Polyethylene Teraphtalate/ Polyethylene (PET/PE) plastic laminate
Product Dimensions:	240 x 135 x 40* mm
Product Weight:	± 38.5 g
Tube Volume:	13 mL
Collected Volume:	8-12 mL
Preservative Volume:	3.4 mL
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Push-fit
Shelf-life:	2 years
Storage temperature:	2-30°C
Sterility:	Non-sterile



\*The height of the pouched device can be squeezed to 30mm when being included in a postal kit