

# Data Sheet

Colli-Pee® containing UCM preservative, FV-5000 series



## Product Description

Product Name:	Colli-Pee® containing UCM preservative, FV-5000 series
Reference Number:	N00176
Components:	Housing, floater, prefilled collector tube, tube cap and IFU
Order Quantity:	25 or 100

## Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of medical grade paper/ multiply PET/PP-plastic laminate
Product Dimensions:	300 x 150 x 46 mm
Product Weight:	± 54.5 g
Tube Volume:	30 mL
Collected Volume:	16.5-23 mL
Preservative Volume:	7 mL
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Screw thread
Shelf-life:	2 years
Sterility:	Non-sterile

