

# Data Sheet

Colli-Pee® containing UCM preservative, FV-5000 series



## Product Description

|                   |  |
|-------------------|--|
| Product Name:     | Colli-Pee® containing UCM preservative, FV-5000 series       |
| Reference Number: | N00176   |
| Components:       | Housing, floater, prefilled collector tube, tube cap and IFU |
| Order Quantity:   | 25 or 100  |

## Product Characteristics

|                                      |   |
|--------------------------------------|---|
| Material:                            | Polypropylene (PP) and High Density Polyethylene (HDPE)                                   |
| Packaging:                           | See-through pouch constructed of medical grade paper/<br>multiply PET/PP-plastic laminate |
| Product Dimensions:                  | 300 x 150 x 46 mm   |
| Product Weight:                      | ± 54.5 g  |
| Tube Volume:                         | 30 mL   |
| Collected Volume:                    | 18-22 mL  |
| Preservative Volume:                 | 7 mL  |
| Collector Tube - Tube Cap Interface: | Screw thread  |
| Housing - Collector Tube Interface:  | Screw thread  |
| Shelf-life:                          | 2 years   |
| Sterility:                           | Non-sterile   |

