

# Data Sheet

Colli-Pee<sup>®</sup>, FV-5000 series



## Product Description

Product Name:	Colli-Pee <sup>®</sup> , FV-5000 series
Reference Number:	N00055
Components:	Housing, floater, collector tube, tube cap and IFU
Order Quantity:	25 or 100

## Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of medical grade paper/ multiply PET/PP-plastic laminate
Tube Volume:	30 mL
Collected Volume:	16.5-23 mL
Preservative Volume:	N/A
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Screw thread
Shelf-life:	5 years
Sterility:	Non-sterile

