

Data sheet

Colli-Pee® FV-5010a-IUO



Product description

Product Name: Colli-Pee® FV-5010a-IUO

Reference Number: N00237_IUO

Components: Housing, floater, cobas PCR Media Tube, tube cap and IFU

Order Quantity: 25 or 100

Product Characteristics

Material: Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging: See-through pouch constructed of Polyethylene Teraphtalate/

Polyethylene (PET/PE) plastic laminate

Product dimensions: 240 x 135 x 40 mm

Product Weight 40.5 g
Tube Volume: 14.5 mL
Collected Volume: 9-12 mL
Preservative Volume: 4.3 mL

Collector Tube - Tube Cap Interface: Screw thread Housing - Collector Tube Interface: Push-fit

Shelf-life: Dependent on shelf-life of cobas PCR media, but max. 5 years

Storage Temperature: 15-30°C Sterility: Non-sterile

