

# Data sheet

Colli-Pee® FV-5004




## Product description

Product Name:	Colli-Pee® FV-5004
Reference Number:	N00290
Components:	Housing, floater, collector tube, tube cap and IFU
Order Quantity:	25 or 100

## Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of Polyethylene Teraphtalate/ Polyethylene (PET/PE) plastic laminate
Product Dimensions:	240 x 135 x 40* mm
Product Weight	± 35 g
Tube Volume:	7 mL
Collected Volume:	3-5 mL
Preservative Volume:	NA
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Screw thread
Shelf-life:	5 years
Storage Temperature:	0-40°C
Sterility:	Non-sterile



\*The height of the pouched device can be squeezed to 30 mm when being included in a postal kit.