

# Data Sheet

## Colli-Pee® Large Volumes



### Product Description

Product Name:	Colli-Pee® Large Volumes
Reference Number:	N00269_RUO
Components:	Housing, floater, collector tube, tube cap and IFU
Order Quantity:	25 or 100

### Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of medical grade paper/multiply PET/PP-plastic laminate
Product Dimensions:	300 x 150 x 46 mm
Product Weight:	± 53 g
Tube Volume:	50 mL
Collected Volume:	33.5-42.5 mL
Preservative Volume:	NA
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Push-fit
Shelf-life:	5 years
Storage temperature:	0-50°C
Sterility:	Non-sterile

