

Data sheet Colli-Pee® FV-5040



Product description

Product Name: Reference Number: Components: Order Quantity:

Product Characteristics

Material: Packaging:

Product Dimensions: Product Weight Tube Volume: Collected Volume: Preservative Volume: Collector Tube - Tube Cap Interface: Housing - Collector Tube Interface: Shelf-life: Storage Temperature: Sterility: Colli-Pee[®] FV-5040 N00269 Housing, floater, collector tube, tube cap and IFU 25 or 100

Polypropylene (PP) and High Density Polyethylene (HDPE) See-through pouch constructed of medical grade paper/multiply Polyethylene Teraphtalate / Polypropylene (PET/PP) plastic laminate 300 x 150 x 46 mm

± 53 g 50 mL 33.5-42.5 mL NA Screw thread Push-fit 5 years 0-50°C Non-sterile

