

Data Sheet

Colli-Pee® Small Volumes



Product Description

Product Name:	Colli-Pee® Small Volumes
Reference Number:	N00290_RUO
Components:	Housing, floater, collector tube, tube cap and IFU
Order Quantity:	25 or 100

Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of Polyethylene Teraphtalate/ Polyethylene (PET/PE) plastic laminate
Product Dimensions:	240 x 135 x 40* mm
Product Weight:	± 35 g
Tube Volume:	7 mL
Collected Volume:	3-5 mL
Preservative Volume:	NA
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Screw thread
Shelf-life:	5 years
Storage temperature:	0-40°C
Sterility:	Non-sterile



*The height of the pouched device can be squeezed to 30mm when being included in a postal kit