

# Data Sheet

## Colli-Pee® Large Volumes containing UAS



### Product Description

Product Name:	Colli-Pee® Large Volumes containing UAS
Reference Number:	N00342_RUO
Components:	Housing, floater, prefilled collector tube, tube cap and IFU
Order Quantity:	25 or 100

### Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of Polyethylene Teraphtalate/ Polyethylene (PET/PE) plastic laminate
Product Dimensions:	300 x 150 x 46 mm
Product Weight:	± 66 g
Tube Volume:	50 mL
Collected Volume:	33.5-42.5 mL
Preservative Volume:	11.5 mL
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Push-fit
Shelf-life:	1 year
Storage temperature:	15-30°C
Sterility:	Non-sterile

