

Data Sheet

Colli-Pee® UCM® FV-5003b



Product Description

Product Name:	Colli-Pee® UCM® FV-5003b
Reference Number:	N00383
Components:	Housing, floater, prefilled collector tube, tube cap and IFU
Order Quantity:	25 or 100

Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of Polyethylene Terephthalate/ Polyethylene (PET/PE) plastic laminate
Product Dimensions:	240 x 135 x 40* mm
Product Weight:	± 26 g
Tube Volume:	3.5 mL
Collected Volume:	1.7-2.9 mL
Preservative Volume:	0.9 mL
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Screw thread
Shelf-life:	9 months
Storage temperature:	2-30°C
Sterility:	Non-sterile

*The height of the pouched device can be squeezed to 30mm when being included in a postal kit

